

## AMORET TRADING 68

LEADING THE WAY IN HEALTHCARE PRODUCTS

## PHYSIOTHERAPY EQUIPMENT

WWW.AMORET.CO.ZA



# CONTACT DETAILS

TEL 013 244 2119 FAX 086 528 7895 CELL 072 446 4043

STEFUS@AMORET.CO.ZA

3 BOSBERG CRESCENT, AERORAND, MIDDELBURG, MPUMALANGA, SOUTH AFRICA

POSTNET SUITE MW177 PRIVATE BAG X1838 MIDDELBURG, 1050

WWW.AMORET.CO.ZA



## **INDEX**



## Page:

Agility Training	1
Bracing Balance	1 1
Exercise and Rehabilitation	2
Gross Motor Development	2
Hand Therapy	3
Hot and Cold Therapy	3-4
Light Therapy	4-5
Lymph Drainage	5
Needling	5
Respiratory	6
Shockwave Therapy	7
Suction Units	7-8
Strapping & Tapping	8
Tens & EMS	8-9
Traction Units	9-10
Treatment Equipment	10-11
Ultrasound / Combination / Lazer Unit	11-12



## **AGILITY TRAINING**



AMORET TRADING 68

LEADING THE WAY IN HEALTHCARE PRODUCTS













A wide variety available on request















## **BALANCE**

A Wide variety available on request











WWW.AMORET.CO.ZA STEFUS@AMORET.CO.ZA

## EXERCISE AND AMORET TRADING 68 REHABILITATION THE WAY IN HEALTHCARE PRODUCTS



## **GROSS MOTOR DEVELOPMENT**




## **HAND THERAPY**



A Wide varie

TheraBand

The





## **HOT & COLD THERAPY**

## **Thermalator (Hydrollocator)**



- Whitehall Thermalators maintain a constant water temperature and come standard with exposed heating elements.
- These units provide for comfortable and effective moist heat therapy, promoting both easy cleaning and safety.
- The tank body, lid, bottom and tank housing are all made from type 304 stainless steel for years oftrouble-free use.
  - Temperature sensor.

### Keno Wax Heater 500/800gm



- This Wax Heater comes with the hot wax insert and a handle for the insert. This unit is suitable for strip wax, hot wax and paraffin wa
- .



## **LIGHT THERAPY**

#### **Photizo Sport**



- Photizo® Sport is effective in treating pain and inflammation associated with overuse injuries, e.g. runner knee, shin splints, tennis and golfers elbow etc.
  - Photizo® Sport stimulates healing, relieves pain, swelling
- Inflammation, reduces muscle spasms, and increases the blood and lymph vessel diameter, to ensure clearing of waste products and inflow of oxygen and nutrients in the areas treated

### **Photizo Pain Relief**



- It is useful for treating pain and inflammation associated with joint diseases such as all types of arthritis and gout.
- Photizo® Pain Relief reduces muscle spasms, and increases the blood and lymph vessel diameter, to ensure clearing of waste products and inflow of oxygen and nutrients in the areas treated.

#### Photizo Light Professional®





- Photizo® Light Professional all in one use.
- Photizo is available for use by professionals such as

Physiotherapists and Veterinarians (Photizo Professional)

• NOW also for the home user in a small, portable unit which is

easy to use (Photizo Home Care).

## LYMPH DRAINAGE

### **Lympha Press Optimul**



- 12 Channel Professional Lymph Drainage System
- · Variety of garments and expanders
- It can run two garments at same tim



## **NEEDLING**

#### SIZES AVAILABLE:

13mm x 0.22mm, 25mm x 0.25mm, 25mm x 0.30mm, 25mm x 0.35mm, 30mm x 0.25mm, 40mm x 0.30mm, 50mm x 0.30mm, 50mm x 0.35mm, 75mm x 0.30mm,



## RESPIRATORY

100mm x 0.30mm, 100mm x 0.35mm;



### **Ultrasonic**

Nebulisation capacity : 3ml/min max

Medicine chamber capacity: 100ml

Particle size : 1 - 5µm

• Dimensions: 30(L) x 15(W) x 21(H)cm

Weight: 2.8kgTimer 0~60min

#### **Hospital Grade Ultrasonic**





Nebulisation capacity : 4ml/min

Particle size : 1 - 5µm

Medicine cup capacity: 50ml max

• Dimensions: 30cm(L) x 18cm(W) x 27cm(H)

Weight : 4.2kgTimer 0-60mins

### **Power Massager Thrive**



This massager combines a vertical motion with an irregular rotary motion which raises massaging effects. This product can clear blood vessels, the twigs of the capillary vessels in particular. Blood circulation will be promoted immediately and vibration relives your muscle from fatigue . The vibration also reduce Phlegm in your lungs. 3 interchangeable heads soothe all areas of the body and vibrating massage. Speed 2500rpm – 3800rpm

## SHOCKWAVE THERAPY



Action Head: Radial

Dimensions: 50x102, 7x46 cm

• Control via touch screen graphics technology 12,1", resolution 600x800 pixels

Patient database with expandable memory

• Encyclopaedia treatments dictionary pathologies, anatomical illustrations and photographic guide to treatment

• Set frequency: 1-22 Hz

• Setting pressure: 1-5 Bar

Maximum positive pressure: 28 MPa

Energy density: 0.68 mJ / mm

• Number of shock therapy: 0-9.999 shocks

• It has two independent outputs allowing connection of two applicators

Includes 4 sets of change heads

#### **Shockmaster 300**



LTHCARE PRODUCTS



Small and lightweight
Silent built-in compressor
Perfect viewing angle for you and your patient
Parameters can be set quickly and easily using the large operating buttons
Clear representation of the treatment time on a colourful illuminated bar
Safe and simple way of putting the hand piece down:
protective mat

## **SUCTION UNITS**

#### Askir 30



- Mobile and portable with a weight of ±3.6kg
- Also available as table device
- Mains-powered device
- Used for minor surgeries as well as nasal, oral and tracheal aspiration
- Maximum 120 minutes operation time, then 60 minutes cooling time
- Ideal for the home as it is light to transport
- Satisfies all current norms and EU safety standards
- Comes with a 1 litre autoclave able plastic collection jar, overflow valve, suction control, vacuum display, connection tube, suction catheter and antibacterial filter
- The stand model comes with 2 x 2 litre autoclave able plasticcollection jar

Hospivac 350



- New Hospivac 350 is a surgical aspirator working with 230V ~ / 50 Hz network electricity
- Four antistatic wheels with two brakes and a movement handle
- Easily moveable and designed for virtually continuous use
- The product is supplied with 2 x 2 litre polycarbonate autoclave able jars with overflow valve
- Equipped with aspiration regulator and vacuum indicator located on the front panel
- Power consumption: 230VA
- Fuse 1 x 4 A 250V
- ISO 10079-1 High Vacuum/High Flow
- Maximum suction aspiration: -0.90 Bar
- Maximum suction flow: 60 l/min
- Dimensions: 460mm(l) x 600mm(h) x 420m



## STRAPPING & TAPING



## **TENS & EMS**

#### **Itouch Easy**



- Seven Preset programmes
- Dual Channel with independent intensity controls
- Powerful output 90mA
- Large LCD Display
- Safety interlock with LCD indicator
- Integral treatment timer and Patient Compliance memory
- Low battery indicator
- Neck cor

### Tenscare T.E.N.S Sport 2



WWW.AMORET.CO.ZA STEFUS@AMORET.CO.ZA

- Li-ion mobile phone style battery
- Power 99mA into 1kOhm enough for the largest muscles
- 2 fully adjustable manual muscle stimulation programmes
- Active rest in muscle stimulation programmes
- 10 preset and 2 manual TENS programmes
- 10 massage programmes
- Special massage programmes for swelling reduction
- Detailed usage memory



#### Socks



 The conductive sock provides stimulation to the entire foot and ankle to increase blood circulation, reduce swelling and relieve pain.

<u>Hand</u>



• The conductive gloves are excellent for use in treatment of Carpal Tunnel Syndrome, arthritic pain and post-operative swelling reduction.

## **TRACTION UNITS**

#### **Vertetrac Portable 3 D Traction**



- Vertetrac's lumber ambulatory decompression system is an innovative and cost-effective approach to treating Lower Back Pain and Sciatica
- Vertetrac aims to bring relief to people suffering from lumbar spine related injuries and disorders, including herniated discs, lower back pain, and Sciatica as an alternative to surgery, body braces, or decompression table treatments
- Outperforms traditional traction by combining decompression with horizontal force to enable simultaneous 3D decompressive traction therapy

## **Mattress Traction Kit**



- Aluminium frame
- Included with 9kg weight bag, 33cm spreader bar, traction cord and traction bracket
  - Adjustable Velcro head halter and a steel pulley support

WWW.AMORET.CO.ZA STEFUS@AMORET.CO.ZA



#### **Overdoor Traction Kit**



 Overhead traction hanger with steel frame. Complete with head halter, spreader bar, cord and water weight bag

## TREATMENT EQUIPMENT

#### **Fusion 2 Section Electric Bed**



- 5cm High density foam
- Height adjustable by Hand or Foot switch
- One Electric Linak Motor
- Retractable Wheels
- Manual backrest +70°
- 3 years warranty for fram

#### **Fusion 5 Section Electric Table**



G2 Duo

- 5cm High density foam
- Height adjustable by Hand or Foot switch
- One Electric Linak Motor
- Retractable Wheels
- Manual backrest +70°
- 3 years warranty for frame



- The basis forcomfortableand safereceipt ofyour clients.
- The longback piecesupports asitting position.
- Yourclient isoptimuminvolved in thetreatment.




## ULTRASOUND / COMBINATION/ LAZER UNITS

#### **Astar Sonaris Ultrasound Unit**



- $\bullet$  S version co-operates with two types of ultrasound heads:
- Effective radiation area of 4 cm<sup>2</sup>, generating ultrasound wave with frequency of 1 MHz or 3,5 MHz,
- Effective radiation area of 1 cm<sup>2</sup>, generating ultrasound wave with frequency of 1 MHz or 3,5 MHz.
  - M version co-operates with one type of ultrasound head:
    - Effective radiation area of 4 cm<sup>2</sup>, generating ultrasound







- 1 and 3 MHz frequencies
- Pulsed and continuous therapy operation (10%, 20%, 50% and 100%)
- Variable duty cycle outputs of 16 Hz, 48 Hz, or 100 Hz
- Clinical indications provide the most appropriate application and parameter set for the patient's condition

### Chattanooga Intelect Mobile Combo



- Variable duty cycle outputs of 16 Hz, 48 Hz, or 100 Hz for ultrasound
  - Multiple waveforms on board
- Clinical Indications on board, providing the most appropriate parameter set for a variety of conditions
  - 15 user-defined memory positions
- Lightweight design, along with battery powered option and custom designed mobile carrying bag, make transport ease

#### **Dundex**





• Combined ultrasound with two electrotherapy channels including interferential traditional (4 pole) and interferential pre modulated (2 pole) functions, which makes the device remarkable in physiotherapy, rehabilitation and adjacent areas

## **Laser Chattanooga Intelect Mobile**



- Class 3B Low Level Laser Therapy
- Independent control over all parameters
- Continuous and pulsed operation available
- Fully variable pulsed frequencies
- Real time feedback of energy delivered
- Display of energy in choice of joules or joules/ cm<sup>2</sup>
- Helpful menu of Clinical Indications specific to the laser applicator



# CONTACT

TEL 013 244 2119 FAX 086 528 7895 CELL 072 446 4043

STEFUS@AMORET.CO.ZA

3 BOSBERG CRESCENT, AERORAND, MIDDELBURG, MPUMALANGA, SOUTH AFRICA

POSTNET SUITE MW177 PRIVATE BAG X1838 MIDDELBURG, 1050

WWW.AMORET.CO.ZA

