

# OCCUPATIONAL THERAPY

WWW.AMORET.CO.ZA

## CONTACT DETAILS



## CHÉRI-LEE BOSHOFF

TEL 013 244 2119 CELL 086 528 7895 FAX2MAIL 086 477 2680

## CHERILEE@AMORET.CO.ZA

3 BOSBERG CRESCENT, AERORAND, MIDDELBURG, MPUMALANGA, SOUTH AFRICA

POSTNET SUITE MW177 PRIVATE BAG X1838 MIDDELBURG, 1050

WWW.AMORET.CO.ZA





## Product Range

3
4
5
6-7
8
9
<u>10</u>
11-14
<u>15</u>
<u>16</u>

## **Strapping Material**





#### Velfoam #2 Soft Padded Loop.

Velfoam has a foam backing to help prevent slippage. Velfoam #2 has hook on two sides making it more versatile than Velfoam. Commonly used for disposable splints. For a stretch-free alternative, try Beta Pile II Loop. Hand wash only. prevents

4.6 m L. Beige 5.1cm



#### Beta Pile II Loop.

Softest padded loop with no harsh surfaces to cause skin irritation. Best for short-term splints (one-two months); upper extremity splints (resting only); paediatric splints and non-Ambulatory AFOs. Stabilizing fabric

deformity. Hand wash only.



Fabrifoam NuStim Wrap

Knitted blend of Nylon and

Lycra with

Patended foam lining.

**Replaces Neoprene wraps** 

for wrist, elbow,

Ankle and knee support.

Also secures proper

placement

Of electrodes.

#### **Co-Flex is Self-Adhesive** No clips or fasteners needed. Maintains Consistent compression with nonslip support. **Contains Latex**



Neoprene Hook Sensitive Fabric on one side, Nylon tricot on the other.

> Standard Hook & Loop Stock Sizes: 25 & 50mm (white) Non-stock sizes: 16, 20, 38 and 100mm Packaging: 25m roll



Self-adhesive Hook & Loop Stock Sizes: 25 & 50mm (white) Non-stock sizes: 16, 20, 38 and 100mm Packaging: 25m roll



#### Webbing

All-purpose webbing provides versatility for splint or equipment modification<sup>4</sup>

WWW STEFUS@AMORET.CO.ZA

## **Padding Material**





Aliplast 4E Medium Soft Grade. Use for splint and cast Padding and resting splint liners. 0.32 cm, Self-Stick, 45.7cm X 60.9 cm

> Alimed Soft Sponge Cushions styloids and bony prominences.



Reduces Shear and abrasion. Use by itself

Or with liner to insulate against

shock and vibration.



#### AliMed<sup>®</sup> Molestick

Molestick is also slightly thicker than Moleskin, for more cushioning. Foam layer compresses underweight to equal Moleskin in thickness. Beige.

> Alimed Moleskin is a soft thick, adhesive Backed cotton cloth that protects exposed skin against painful blisters, calluses and spots.





Plastazote Self Adhesive Padding is made from

Plastazote. Soft Feel, can be applied to splinting

#### Neoloop<sup>™</sup> Material

Ideal for fabricating stiffer custom positioning splints, devices or padded handles.

WWW.AMORET.CO.ZA STEFUS@AMORET.CO.ZA

soft material before forming

## Scar Tissue Management

### **Elastomers**

Help control and manage hypertrophic scarring. Use elastomer to form total-contact skin models to be worn with burn garments Gloves and hand splints. Use with stockinette to form soft splints

#### Pediplast includes catalyst

#### Otoform







#### Gel Dots and Squares 3-mm thick

polymer gel for versatile protection against pressure and friction.





**Crystal™ Clear Polymer Gel** Crystal Clear Polymer Gel contains pressure-sensitive, medical-grade mineral oil adhesive backing. Apply gel to skin. Pressure- and shear-absorbing interface. Cuts to size. Contours. 1/8" x 4" x 4". (3.2mm x 10.2cm x 10.2cm)



**Pressure Garment Material** Elastonet - Ealstic Spandex/Nylon material for pressure garments for the treatment of burns, scars and edema control.



#### **Gentle Compression Glove**

Offers outstanding value for edema control. Gloves offer equivalent performance at a reasonable price.



#### Mini massager

Ideal for scar massage, desensitization, muscle stimulation, oral stimulation and sensory

### <u>Ortho</u>





#### DeRoyal® LMB Spring Finger Extension Assist (Model 501

(Model 501) extends PIP joint with slight extension on MP. Adjustable. For finger flexion tightness and boutonniere (allows IP flexion). Force at 15° flexion is approximately 8 oz.; 90° flexion is approximately 4 lbs. Wire-Foam<sup>™</sup> for a customized fit

Bunnel/Reverse Knuckle Bender Extends PIP joint. Black felt with Latex-free rubber bands.





#### Alimed Dynamic Digit Extensor Tube It is the better way to reduce PIP and

DIP joint finger flexion contractures. A curved banana shaped, neoprene tube locates a 30 degree reverse angle under the PIP joint, stimulating a spring-like mechanism.

#### Alimed Antimicrobial Blue Carrot

It has a soft covering permanently impregnated with Agion Antimicrobial odour control. Blue version solves the toughest hand Contracture problems.



#### **Oval-8® Finger Splints**

Oval-8 Finger Splint - shown for mallet fingerOval-8 Finger Splints are the ideal choice whenever you need to stabilize or immobilize the small joints of your patient's fingers.

#### D-Ring Regular Thumb Splint

Protects and supports the thumb and wrist for conditions such as arthritis, tendinitis, sprains, strains, skier's thumb, gamekeeper's thumb, de Quervain's syndrome, wrist instabilities and postsurgery





#### Alimed Ball Splint

Uses a reflex inhibiting position to control spasticity and tone. Adjustable thermoplastic.



#### **Copper gloves**

Instant relief: Compression gloves made with real copper reduces inflammation and swelling, making them perfect arthritis gloves.

#### Freedom Swedish AFO

For Static dorsiflexion assistance and lateral stability for entire foot/ankle area. Injection moulding allows for thicker polyethylene on vertical aspect for rigidity

#### Silipos Digistrip

A half-coated gel lining that provides more extensive gel coverage Around the affected area to help relieve pressure and friction. It is Best used on blisters, corns, hammertoes, calluses and nailbed injuries.







#### Viscolas Solid Polyurethane Heel Spur Cushion

Viscolas is a solid urethane polymer insoling material precisely Tailored for heat, high stress and repetitive impact unlike Sorbothane. Viscolas insoles have an extremely long life.



## Splinting Material MAIDFIRM



M sheets (transparent when heating)

E sheets( not transparent when heating)


	M Thermoplastic sheet (P sheet)				E Thermoplastic sheet (K sheet)						
Thickness	1.6mm	2.0mm	2.4mm	3.2mm	4.0mm	3.2mm, 4.0mm					
Perforation	0%	1%	5%	12%		1%	5%	12%	22%		
Feature	Transparent, 100% memory,good tensile, light weight					-	high strength,good tensile,light weight				
	White, Beige,Blue, Red Special color, Non- perforated, Half piece/special size are available										

Our low temperature thermoplastic is made of eco-friendly polycaprolactone materials from UK, it's strong, rigid immobilization and easy moldability(activated at 65°C-70°C), with less shrinkage, with zero tackiness, no stick to skin, hair.

#### CE, FDA & ISO 13485 approval









**Klarity White** 

Klarity White is similar to Aquaplast and has 100% memory turning transparent when completely heated.



#### Klarity KS

Klarity KS has maximum rigidity and is ideal for fracture bracing and trunk immobilization.



#### **Contour Blend**

Highly conforming and fast setting. Contour blend is our most versatile material.



#### **Klarity Easy Fit**

Klarity Easy Fit makes strong rigid splints for patients who are hard to position or have strong joint contractures.



#### **Contour colours**

are ideal for circumferential splints, neurological splints, all hand and wrist splinting.



#### **Contour Pro**

Contour Pro makes the most rigid and precise fitting splint and orthosis.

STEFUS@AMORET.CO.ZA

www

## **Splinting Accessories**



#### **Resolving Hole Punch**

Self-Opening, Nickel-plated, Punches from 0.2 cm to 0.48cm Throat 8.3 cm



#### **Beta Pile II Self Adhesive Straps**

Combine stretch-resistant padded Beta Pile With 7.6cm of self-adhesive hook-and-loop. Trim loop to length. Won't stretch. Beige

**Double Sided Hooks** 



**Zytel D-Rings** for fabrication Custom D-Ring straps. 2.5 cm. 100/pk



**Fiskars Scissors** 

**Splinting Scissors** 



Beta Pile II



D-Ring Straps combine soft Beta Pile II With 5.1 cm of hook-and-loop closure. Nylon D-Ring. Comfortable for short-term Wear. Cannot be trimmed.




#### Halfsize Thermal Bath

Dual Thermostat is factory set at C (High) and C low. Hinged insulated covers. Stainless Steel. Rear Drain. Finger Loops Have a metal eyelet. 10/pk

**Finger Hooks** 

**Speed Rivets** 



## Hand Exercise Equipment




Alimed Color-Coded Therapy Putty It comes in resistances consistent with industry standards. Silicone based putty





Economy Hand Exerciser Add bands to increase tension.





#### Alimed T-Foam Cubes

For use in early joint mobilization exercises, these slow-recovery cubes help patient safely pace prescribed exercise activity. It returns to its original shape when released.



#### **Digi-Flex**

Develops isolated finger strength, flexibility and coordination as it builds hands and forearm strength.



WWW.AMORET.CO.ZA STEFUS@AMORET.CO.ZA Maintain ROM and build upper extremity strength in the hemiplegic arm. Patient follows figure 8 pattern to complete a sequence of motions. Use without resistance or add weight.

## Aids for Daily Living





Folding Commode Grey powder-coated steel. Plastic Arms, commode bucket, cover and splash shield included.



Ambulation Belt Is wider than typical gait belts. Nylon web

handles and metal



#### Bathtub Rail

Fits most tubs. Tool-free installation. Vinyl-coated, heavy gauge welded steel. Ring rubber pads will not mar.





WWW.AMORET.CO.ZA STEFUS@AMORET.CO.ZA

#### Patient Handler Strong. Flexible poly-urethane ensures safe and secure patient transfer that is less stressful on the attendant's back.

Safety NoSlip Mat

Only 3.2cm thick so wheelchairs can roll right over it. Offers traction wherever there are slippery surfaces.



Dycem is a no slip Rubber-like plastic, material used to stabilize. Reusable. Cuts with scissors. Cleans with soap and water. Blue. Matting is 1/32" thick Pads are 3/16" thick Not made of natural rubber latex. Long lasting.



#### Raised Toilet Seat with Armrest

It features a twist-and-holding mechanism. All plastic seat and clamp will not corrode. Padded armrests assist sitting and standing.



#### Tub/Bath Grab Bars Include everything needed to mount on any surface.





**Aluminium Four-Wheel Roller** 

Frame: Aluminium, 5896.7g, 136,078g capacity Locking Loop Breaks, 15.2 cm Wheels, Accessories Basket/Back Rest 83.8 cm-96.5 cm Handle Height Adjustment

> Economy Molded Tray 100% Recycled Plastic. Various Colours



Heavy Duty Wheelchair Gloves Foam padded black leather.





#### Lift Away Tray

Secures to the wheelchair armrests with hook and loop strips Tray helps achieve proper body positioning and alignment. Comfortable foam padding. Wipe-clean vinyl.



Alimed Plus Folding

Same feature as the

but it folds for easy

storage or travel.

Alimed Plus Reacher,

Reacher





**Hip Kit** The ADL items most often recommended by therapists for post-hip surgery patients. 66 cm UltraLight reacher,

Molded sock aid with two cords and foam handles



Makes bathing and dressing easier for anyone with limited reach or mobility



**Economy Hip/Knee Kit** For the patient recovering from total hip or knee replacement surgery.



No-Tie Stretch Lace with pull-tab lock.



**Hook-End Plastic** Shoehorn



**Ableware Button** Hook/Zipper Pull Soft, cushioned handle for comfortable grasp



#### Dress EZ

With large 2.5 cm-diameter soft foam handle. Serves as a dressing aid and a shoehorn



**Universal Elastic ADL** Cuff slides on for a secure, strapless fit. Slide eating utensils into fixed pocket to help residents eat without assistance.

### Grip Utensils

These easy togrip, high quality utensils have comfortable, oversized, textured handles. Stainless steel.



3 68



Colour-Coded Handle Foam Is slipresistant, nonabsorbent and dishwasher safe.

Rocking T Knife Requires less strength and dexterity than other rocker knives.

#### Flexi-Cut Cups Generously cutout rims. Control and direct liquid flow by squeezing cups. Dishwasher safe.



Large Plastic Plate Guard

Dishwasher

safe.



**STA-Pen Writer** Glides along the paper to facilitate writing for people with limited hand strength or control



Weighted Universal Holder Slips onto pens, pencils, markers, paint brushes and more.



**ETAC Food Prep Board** Helps user open a jar, grip a bowl, peel, or cut with one hand. Quick clamp holding.


#### WWW.AMORET.CO.ZA STEFUS@AMORET.CO.ZA

16



## Assessment Tests

(Wide range available please contact us for more information)



**Heydes Modapts** 




DLOTCA



Cognitive Assessment of Minnesota





Developmental Test of Visual Perception



DTVP-A



Rivermead Behavioural Memory Test-Third Edition





**Spinal Function Sort** 



Hand Function Sort

## Assessment Tools



Baseline Complete Hand Evaluation Set All the instruments you need for a complete hand evaluation stored in a portable carrying case.



#### **Personal Rulangemeter** Has a 360° dial in 5°

increments. Plastic with linear measurement to 15.3cm

**Personal Goniometer** Has two 0° to 180° (5° increments) scales that read in opposite directions.





#### Dynamometers

Frequently used grip strength measurement instruments with heavy-duty handles and five grip positions. Easy to adjust. Record in both pounds and kilograms.



#### Perdue Pegboard

Includes instructions and standards for evaluation, pegboard, pegs, collars and washers O'Conner Tweezer Dexterity Test Includes board, pins and manual with norms.



## AMORET TRADING 68



Baseline LTE Three Piece Hand Evaluation Set Includes all instruments necessary to evaluate strength and range of motion of hand and fingers in a portable carrying case.



Baseline Stainless Steel Goniometer

**O'Conner Finger Dexterity Test** Measures patients small-part manipulation ability.



WWW.AMORET.CO.ZA















## **Suppliers**

We choose to only represent international companies which are on the forefront of the medical profession.

AliMed, INC.



A leading manufacturer of and supplier of high quality products - everything from preformed orthoses, alarms for fall management, cushions, diagnostic imaging, and operating room accessories, orthopaedic rehabilitation equipment, and ergonomic workplace solutions.



>>> For more information: www.alimed.com

Maidfirm medical devices Co Ltd.



Manufacturer of low-temperature thermoplastic splinting material, thermoplastic radiotherapy positioning fixed materials, rehabilitation and orthopaedic external fixation materials for occupational therapy and orthotist applications.

>>>> For more information: <a href="http://www.maidfirm.com">www.maidfirm.com</a>

#### **Klarity Medical Products INC**

Klarity

Manufacturer of low-temperature thermoplastic splinting material for occupational therapy and orthotist applications (USA).

>>> For more information: www.klaritymedical.com

#### 3-Point Products™

Innovations by design! A manufacturer of innovative and comfortable, attractive, and affordable hand, finger wrist and elbow splints.

>>> For more information: <a href="http://www.3pointproducts.com">www.3pointproducts.com</a>

**VALPAR International Corporation** 



VALPAR International manufactures and develops specialized test instruments and software that measure workrelated skills. Their products are used worldwide in education, workforce development, allied health, and are the industry standard for work skills assessment.

>>> For more information: <a href="http://www.valparint.com">www.valparint.com</a>



#### **BASES of Virginia**



Busic and Scott Evaluation Systems of Virginia (BASES of VA) came into being in late 2008 after VALPAR Int sold its rights for the manufacturing of some of the Component Work Samples. They are an authorized manufacturer & distributor of Valpar Component Work Samples (VCWS). VCWS have earned an excellent reputation since their initial introduction in 1973.

Work Samples are generalized work-like tasks that are administered under specific instructions. The challenge and gamesmanship that the Work Sample presents stimulates the client to put forth their best effort. Language and reading requirements are low, so they do not present barriers. The scores are interpreted in a variety of ways, including MTM standards, DOTWorker Qualification Profiles and norms.

>>> For more information: <a href="http://www.basesofva.com">www.basesofva.com</a>

## CHÉRI-LEE BOSHOFF

TEL 013 244 2119 CELL 086 528 7895 FAX2MAIL 086 477 2680

## CHERILEE@AMORET.CO.ZA

3 BOSBERG CRESCENT, AERORAND, MIDDELBURG, MPUMALANGA, SOUTH AFRICA

POSTNET SUITE MW177 PRIVATE BAG X1838 MIDDELBURG, 1050

WWW.AMORET.CO.ZA





LEADING THE WAY IN HEALTHCARE PRODUCTS